




Model # 1200-R Stage

Shown with Optional Rotator #ROT-1200



- **16 lbs.**
- **Rack and Pinion Drive**
- **Smooth Ball Bearing Travel**
- **Custom Holes Available**

- **Black Anodized**
- **Stainless Steel Raceways**
- **One Year Warranty**
- **Co-Axial Knobs**



PROPRIETARY
 This document contains information that is the sole property of Franklin Mechanical & Control Inc. and cannot be copied, reproduced or utilized in any way without the express written consent of Franklin Mechanical & Control Inc.

Franklin Mechanical & Control, Inc.
 8901-B Murray Ave., Gilroy, CA 95020 408-842-0313 www.franklinstages.com