

INFINITY CAPTURE is an intuitive user interface which includes all of the basic features needed to control your INFINITY camera and capture images.



INFINITY CAPTURE

IMAGING SOFTWARE




Lumenera's entry-level software is included with all INFINITY USB 2.0 cameras.

Features include:

- Real Time Video Preview
- Separate Preview and Capture Controls
- Automatic/Manual White Balance
- Manual Hue and Saturation
- Manual Preview Quality
- Multiple Image Averaging
- Global Electronic Gain
- Manual Contrast/Brightness/Gamma
- Automatic/Manual Exposure
- Overlay Bitmaps
- Flat Field Correction
- Highlight Saturated Pixels
- Manual Preview Zoom Level
- Preview Histogram
- User Defined Light Source
- Alpha Blended Bitmaps

Additionally, a USB 2.0 plug-and-play interface removes the need for framegrabbers and allows fast video. Included is a TWAIN interface, allowing the camera to be used with 3rd party software applications.



WWW.OEM-OPTICAL.COM
WWW.OEM-OPTICAL.NET