TLC Techni-Light®



Features

- Provides Even & Diffused Lighting
- Hi-Frequency Output
- · Accurate Color Rendition
- Durable
- High Quality
 Fluorescent Lighting

TIC

TechniQuip offers three new durable low cost Techni-Lights® to assist your microscope illumination needs (TLC -3, TLC-4, and TLC - 5). Each of our Techni-Lights® is a perfect choice for applications requiring a high quality fluorescent illuminator at an economical price. Each illuminator is ideal for microscopic and macroscopic inspection, quality control, photographic applications, assembly and biological analysis. Techni-Light® bathes the subject with 360 degrees of shadow- free illumination. Our reflector housing is cast from high quality Al. alloy for durable and cool operation. Techni-Lights are treated with a baked on enamel which will provide long-term durability.

Manufacturers of Quality Industrial and Medical Instruments



Optional ESD (Electronic Static Discharge) Coating.

ESD Option Is Available In Black Only.

Detachable Modular Cable

Instant Start Electronic Ballast

TLC-4 Techni-Light* Features

Instant Start Electronic Ballast

High Frequency

Variable Output

Intensity Adjustments

Regulated Power Unit

Detachable Modular Cable (3" standard- longer length available)

TLC-5 Techni-Light* Features

Optional ESD (Electric Static Discharge) Coating. ESD Option Is Available In Black Only

Fixed Output

Plugs Into Wall (Optional 240 Vac Ballast is "in line")

The Techni-Light®, exclusively from TechniQuip, is the system for providing even diffused illumination . With the WTL lamp the color temperature is 6000°K simulating natural white daylight. Special electronic circuitry makes TechniQuip's Techni-Light's ideal for frame grabbers and machine vision system.

Each Techni-Light® is treated with a baked on enamel which seals against particulate contamination making it useful in clean rooms. Techni-Lights® come with rubber tipped thumbscrews that provide simple and versatile mounting eliminating the need for adapters in most situations. These Techni-Lights clamp on to lenses from sizes 1.25" to 2.15", and up to 2.6" with set screws. Special and custom threaded microscope adapters are available. Each Techni-Light® has a modular construction which permits easy repairs and lamp replacement. All Techni-Light models are backed by a 2 year parts and labor warranty (excluding lamp). Techni-Light's* lamp life is over 4600 hours.



Special and OEM

configurations

readily available.







